

RayBio® Human Angiogenesis Antibody Array 1

Detect 20 Angiogenic Factors in One Experiment

	A	B	C	D	E	F	G	H
1	POS	POS	NEG	NEG	Angiogenin	EGF	ENA-78	b FGF
2	POS	POS	NEG	NEG	Angiogenin	EGF	ENA-78	b FGF
3	GRO	IFN- γ	IGF-I	IL-6	IL-8	LEPTIN	MCP-1	PDGF-BB
4	GRO	IFN- γ	IGF-I	IL-6	IL-8	LEPTIN	MCP-1	PDGF-BB
5	PIGF	RANTES	TGF- β 1	TIMP-1	TIMP-2	Thrombopoietin	VEGF	VEGF-D
6	PIGF	RANTES	TGF- β 1	TIMP-1	TIMP-2	Thrombopoietin	VEGF	VEGF-D
7	BLANK	BLANK	BLANK	BLANK	BLANK	BLANK	Neg	POS
8	BLANK	BLANK	BLANK	BLANK	BLANK	BLANK	Neg	POS

* For use with serum, plasma, condition medium, urine and all other body fluids, cell lysates and certain tissue lysates samples