



RayBiotech, Inc.

3607 Parkway Lane suite 200
Norcross, GA 30092
Tel: 770-729-2992, 1-888-494-8555
Fax: 770-206-2393
Website: www.raybiotech.com
Email: info@raybiotech.com

Certificate of Analysis and DataSheet Goat Anti-8-Hydroxyguanosine

Catalog No.
MD-19-0003

Isotype
N/A

Description:	Goat A' 8-Hydroxyguanosine Goat Antibody to 8-Hydroxyguanosine
Specificity:	8-Hydroxy-2'-deoxyguanosine and 8-hydroxyguanosine
Host Animal:	Goat
Immunogen:	8-Hydroxyguanosine-conjugate
Format:	Neat, Liquid
Purification:	Not applicable
Concentration:	Not determined
Buffer:	Not applicable
Preservative:	0.1% Sodium azide
Applications:	Suitable for use in ELISA using immunizing peptide (>1:250,000) and IHC using formalin-fixed, paraffin embedded sections (1:200). Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Storage:	Short-term store at 2-8°C. Long term store at -20°C. Avoid multiple/freeze thaw cycles.

**The products are furnished for LABORATORY RESEARCH USE ONLY.
Not for diagnostic or therapeutic use.**