



## Adipokine Related Recombinant Proteins

<u>SPECIFICITY</u>	<u>HOST</u>	<u>DESCRIPTION</u>	<u>RayBiotech Catalog</u>	<u>UNIT</u>	<u>RayBiotech Price</u>
Acrp30	Mouse	Recombinant Mouse Adiponectin glycosilated, HEK cell line	IP-01-435P	2µg	98.00
Acrp30	Human	Recombinant Human Adiponectin glycosilated, HEK cell line	IP-01-434P	2µg	98.00
Acrp30	Human	Recombinant Human Adiponectin His	IP-01-433P	10µg	98.00
Acrp30	Mouse	Recombinant Mouse Adiponectin	IP-01-432P	20µg	98.00
Acrp30	Human	Recombinant Human Adiponectin	IP-01-280P	5µg	98.00
Acrp30	Human	Recombinant Human Adiponectin Globular	IP-01-277P	2µg	135.00
Acrp30	Mouse	Recombinant Mouse Adiponectin Trimeric Form	IP-01-247P	2µg	98.00
Acrp30	Human	Recombinant Human Adiponectin Trimeric Form	IP-01-233P	2µg	98.00
ADFP	Human	Recombinant Human Adipose Differentiation-Related Protein	IP-03-405P	2µg	\$ 98.00
GHRL	Human	Recombinant ghrelin/obestatin preprohormone	RB-08-0001P	10ug	149
GHRL	Human	Recombinant ghrelin/obestatin preprohormone	RB-08-0001P	25ug	253
GHRL	Human	Recombinant ghrelin/obestatin preprohormone	RB-08-0001P	50ug	379
GHRL	Human	Recombinant ghrelin/obestatin preprohormone	RB-08-0001P	100ug	531
Lectin		Lectin	IP-03-545P	2µg	98.00
LECTIN	Plants	NATIVE SOYABEAN SOYABEAN LECTIN:FITC	DS-01-0279	5mg	339
Leptin	Chicken	Recombinant Chicken Leptin Binding Domain	IP-01-509P	10µg	98.00

Leptin	Human	Recombinant Human Leptin Binding Domain	IP-01-508P	10µg	98.00
Leptin	rabbit	Recombinant Rabbit Leptin	IP-01-507P	100µg	378.00
Leptin	dog	Recombinant Dog Leptin	IP-01-506P	100µg	378.00
Leptin	Chicken	Recombinant Chicken Leptin	IP-01-505P	100µg	378.00
Leptin	horse	Recombinant Horse Leptin	IP-01-504P	100µg	378.00
leptin	bovine	Recombinant Bovine Leptin	IP-01-502P	100µg	378.00
Leptin	Mouse	Recombinant Murine Leptin	IP-01-351P	5mg	594.00
Leptin	Mouse	Recombinant Murine Leptin	IP-01-351P	1mg	189.00
Leptin	Mouse	Recombinant Murine Leptin	IP-01-351P	200µg	98.00
Leptin	Ovine	Ovine Leptin	IP-01-239P	100µg	378.00
Leptin	Ovine	Ovine Leptin	IP-01-239P	200µg	475.00
Leptin	Ovine	Ovine Leptin	IP-01-239P	1mg	1445.00
Leptin	Human	Human Leptin	IP-01-228P	1mg	189.00
Leptin	Human	Human Leptin	IP-01-228P	200µg	98.00
Leptin	Human	Human Leptin	IP-01-228P	5mg	594.00
Leptin	Rat	Rat Leptin	IP-01-227P	5mg	594.00
Leptin	Rat	Rat Leptin	IP-01-227P	1mg	189.00
Leptin	Rat	Rat Leptin	IP-01-227P	200µg	98.00
Leptin-qA	Ovine	Recombinant Ovine Leptin Quadruple Antagonist	IP-01-357P	100µg	189.00
Leptin-qA	Ovine	Recombinant Ovine Leptin Quadruple Antagonist	IP-01-357P	20µg	98.00
Leptin-qA	Ovine	Recombinant Ovine Leptin Quadruple Antagonist	IP-01-357P	1mg	1,130.50

Leptin-qA	Human	Recombinant Human Leptin Quadruple Antagonist	IP-01-353P	1mg	1,130.50
Leptin-qA	Human	Recombinant Human Leptin Quadruple Antagonist	IP-01-353P	100µg	189.00
Leptin-qA	Human	Recombinant Human Leptin Quadruple Antagonist	IP-01-353P	20µg	98.00
Leptin-tA	Ovine	Recombinant Ovine Leptin Triple Antagonist	IP-01-356P	20µg	98.00
Leptin-tA	Ovine	Recombinant Ovine Leptin Triple Antagonist	IP-01-356P	100µg	189.00
Leptin-tA	Ovine	Recombinant Ovine Leptin Triple Antagonist	IP-01-356P	1mg	1,130.50
Leptin-tA	Rat	Recombinant Rat Leptin Triple Antagonist	IP-01-355P	1mg	1,130.50
Leptin-tA	Rat	Recombinant Rat Leptin Triple Antagonist	IP-01-355P	100µg	189.00
Leptin-tA	Rat	Recombinant Rat Leptin Triple Antagonist	IP-01-355P	20µg	98.00
Leptin-tA	Mouse	Recombinant Murine Leptin Triple Antagonist	IP-01-354P	20µg	98.00
Leptin-tA	Mouse	Recombinant Murine Leptin Triple Antagonist	IP-01-354P	100µg	189.00
Leptin-tA	Mouse	Recombinant Murine Leptin Triple Antagonist	IP-01-354P	1mg	1,130.50
Leptin-tA	Human	Recombinant Human Leptin Triple Antagonist	IP-01-352P	20µg	98.00
Leptin-tA	Human	Recombinant Human Leptin Triple Antagonist	IP-01-352P	100µg	189.00
Leptin-tA	Human	Recombinant Human Leptin Triple Antagonist	IP-01-352P	1mg	1,130.50
RBP4	Human	as Western blotting positive ctrl**	RB-08-0016-WBC	100ul	69
RBP4	Human	Recombinant retinol binding protein 4, plasma	RB-08-0016P	10ug	66
RBP4	Human	Recombinant retinol binding protein 4, plasma	RB-08-0016P	50ug	199
RBP4	Human	Recombinant retinol binding protein 4, plasma	RB-08-0016P	100ug	299
RBP4	Human	Recombinant retinol binding protein 4, plasma	RB-08-0016P	500ug	448
Resistin	rat	Recombinant Rat Resistin	IP-01-458P	5µg	98.00

Resistin	Mouse	Recombinant Mouse Resistin	IP-01-457P	5µg	98.00
Resistin	Human	Recombinant Human Resistin	IP-01-456P	5µg	98.00
Visfatin	Mouse	Recombinant Mouse Visfatin	IP-08-447P	5µg	98
Visfatin	Mouse	Recombinant Mouse Visfatin	IP-08-447P	25µg	189
Visfatin	Mouse	Recombinant Mouse Visfatin	IP-08-447P	1mg	inquiry
Visfatin	Human	Recombinant Human Visfatin	IP-08-318P	5µg	98
Visfatin	Human	Recombinant Human Visfatin	IP-08-318P	25µg	189
Visfatin	Human	Recombinant Human Visfatin	IP-08-318P	1mg	inquiry
TSH	Human	Recombinant Human Thyroid Stimulating Hormone	IP-02-252P	10µg	\$ 262.50
TIMP3	Human	Metalloproteinases 3	RB-15-0001P	50ug	189