



RayBiotech, Inc.

3607 Parkway Lane suite 200
Norcross, GA 30092
Tel: 770-729-2992, 1-888-494-8555
Fax: 770-206-2393
Website: www.raybiotech.com
Email: info@raybiotech.com

Certificate of Analysis and Data Sheet Goat Anti Carcinoembryonic Antigen (CEA)

Catalog No.
MD-16-0012

Species
Human

Preparation

Host Animal: Goat
Immunogen: Purified Human CEA
Purification: Not applicable

Formulation

Format: Neat, Liquid
Concentration: Not determined
Titer: 1:200,000 (ELISA)
Buffer: Not applicable
Preservative: 0.1% Sodium azide

Storage

Short-term stored at 2-8° C. Long term, aliquot and stored at -20° C.
Avoid multiple freeze/thaw cycles.
Centrifuge before opening to ensure complete recovery of vial contents.

Applications

Table Summary of antibody applications and working conditions

Options / Functions	YES	NO	Not determined	Recommended Work dilution or concentration
ELISA	.			

Note: Other applications are not tested yet. Optimal dilutions should be determined by each laboratory for each application.

**The products are furnished for LABORATORY RESEARCH USE ONLY.
Not for diagnostic or therapeutic use.**