



RayBiotech, Inc.

3607 Parkway Lane suite 200
Norcross, GA 30092
Tel: 770-729-2992, 1-888-494-8555
Fax: 770-206-2393
Website: www.raybiotech.com
Email: info@raybiotech.com

Certificate of Analysis and Data Sheet

Rabbit Anti Human Carcinoembryonic Antigen (CEA)

Catalog No.
MD-16-0015

Species
Human

Preparation

Host Animal: Rabbit
Immunogen: Purified human carcinoembryonic antigen (CEA)
Purification: Not applicable.

Specificity

Human carcinoembryonic antigen (CEA)

Formulation

Format: Neat, Lyophilized.
Reconstitute in 20ul double distilled water.
Concentration: Not applicable.
Buffer: Not applicable.
Preservative: None

Storage

Store lyophilized product at 2-8°C. After reconstitution, store at -20°C.
Avoid multiple freeze/thaw cycles.
Centrifuge before opening to ensure complete recovery of vial contents.

**The products are furnished for LABORATORY RESEARCH USE ONLY.
Not for diagnostic or therapeutic use.**



RayBiotech, Inc.

3607 Parkway Lane suite 200
Norcross, GA 30092
Tel: 770-729-2992, 1-888-494-8555
Fax: 770-206-2393
Website: www.raybiotech.com
Email: info@raybiotech.com

Applications

Table Summary of antibody applications and working conditions

Options Functions	YES	NO	Not determined	Recommended Work dilution or concentration
ELISA			.	
Western Blotting			.	
Immunohistology			.	
EIA			.	
Flow Cytometry			.	
Immunofluorescence staining			.	
RIA	.			1:1,000

Note: Other applications are not tested yet. Optimal dilutions should be determined by each laboratory for each application.

**The products are furnished for LABORATORY RESEARCH USE ONLY.
Not for diagnostic or therapeutic use.**