



RayBiotech, Inc.

3607 Parkway Lane suite 200
Norcross, GA 30092
Tel: 770-729-2992, 1-888-494-8555
Fax: 770-206-2393
Website: www.raybiotech.com
Email: info@raybiotech.com

Certificate of Analysis and Data Sheet

Human Fibrinogen

Catalog No.
MD-20-0060P

Source
Human Plasma

Description:

Purified Human Fibrinogen. Molecular weight: 340 kDa. Isoelectric point: 5.1–6.3

Source:

Human Plasma

Format:

Purified, Liquid

Purification:

Prepared by conventional and affinity techniques

Concentration:

5.3 mg/ml (OD280nm, $E^{1\%} = 15.1$)

Buffer:

10 mM Sodium citrate, 10 mM Sodium phosphate, pH 7.3

Preservative:

None

Applications:

Specific methodologies have not been tested using this product.

Storage:

Store at -80°C . Avoid multiple freeze/thaw cycles.

Inactivation:

Not applicable

**The products are furnished for LABORATORY RESEARCH USE ONLY.
Not for diagnostic or therapeutic use.**